







AccessRx Secure*

Helping to Close the Loop on Medication Administration.

The AccessRx Series of medication delivery carts provides solutions to more securely and efficiently transport and administer medications at the patients' bedside.





	AccessRx MD⁻	AccessRx Secure
Maximum Number of Bins	12	12
Configurable Drawer Options	\checkmark	✓
Individual Locking Patient Drawers	√	✓
Software links to the HIS	✓	
Controlled Substance Drawer Option	✓	✓
Patient Bin Assignment	Automatic	Manual
Antimicrobial Infused Cart Material	MICROBAN	MICROBAN

KEY BENEFITS



Individual Locking Bins & Drawers



Supports real-time BCMA



Supports Centralized or Decentralized Medication Fill Workflow



Exchangeable Medication Cassettes



Secure & Track Patient Medication



Advanced Li-Nano Technology w/ 5 Year Warranty Closing the loop of medication delivery between pharmacy and the patient bedside.

AccessRx MD™



Medication Security & Visibility at the bedside

- Individually locking, patient-specific drawers can be accessed by scanning the patient ID barcode
- Real-time bedside medication administration, barcode patient identification and direct eMAR API access streamlines BCMA process, reduces redundancy and increases accuracy
- Patient drawers automatically assigned using the ADT (Admission/Discharge/ Transfer) interface, eliminating time consuming and error-prone manual drawer labeling
- Utilizes the hospital's Active Directory for consistent user authentication
- Configurable drawer display fields allow customization of the user interface

More Control & Flexibility for Pharmacists across a variety of pharmacy distribution models

- Browser-based AccessManager[™] console provides access to reports while allowing for Active Directory-compatible user management, fleet management and system configuration
- Configurable, time-based access controls reduce the risk of medication administration at an incorrect time and aid in pharmacy reconciliation for unused medications
- Optional medications orders interface Rx MD can act as a patient medication order filler and verify that an order has been administered.
- Secure Pharmacy email alerts for transferred & discharged patients



Streamline Your Workflow

- Secures and tracks patient specific medications for up to 12 patients in individually-locking drawers, reducing trips and time spent at the ADC
- ADT and patient orders transmitted via wireless network to the cart to facilitate faster medication delivery
- Reduce trips to and time spent at the automated dispensing cabinet

Scalable & Centralized Software Architecture and Integration

- Multi-Org scaleable architecture ability to scale from two bed clinic to multi facility health system
- Ability to set up health system access rights, such as by department
- Enterprise reporting available
- Enterprise multi-facility deployment
- HL7 compliant, ADT interface (patient demographics), and RDE interface (Pharmacy approved orders)
- · Configurable interface engine easily communicates with the HIS
- Project management, installation and training by TouchPoint Medical's professional services team

Store and secure your medications to the bedside.

AccessRx Secure









Secure your medications throughout your workflow

- · Individual locking drawers
- · Easy-to-use software
- Storage configurations to secure a variety of medications and controlled substances
- Smoked polycarbonate 4" & 6" drawers
- Aluminum locking storage and controlled substance drawers
- Configurable drawer display fields allow customization of the user interface
- Add AccessRx MD[™] software options to existing AccessRx Secure[™] carts for a scalable, custom solution

152mm (6in)

Configure your system to meet your medication and supply requirements.





DimensionsWidth
Height
Depth
Volume

Drawer Specifications	Drawer Specifications
91.4mm (3.6in)	142.2mm (5.6in)
238.8mm (9.4in)	238.8mm (9.4in)
55.9mm (2.2in)	55.9mm (2.2in)
1213cm³ (74in³)	1885cm³ (115in³)

102mm (4in)

Full Width Drawer Specifications

271.8mm (10.7in)
228.6mm (9in)

50.8mm (2in)

3195cm3 (195in3)

Controlled Substance



Contact your TouchPoint Medical representative to learn more and discuss additional accessories specific to your point-of-care application. Visit **touchpointmed.com/contact** to get in touch with your local representative.

20211027-BR-MedicationDeliveryCarts-EN