



AccessMount™

Reduce Downtime & Extend the Life of Your Equipment

Your Touchpoint Medical service technician can design and implement a maintenance or refresh program to lower your overall support costs, maximize your investment, and keep your teams focused on patient care.

Preventive maintenance can include replacement of worn bearings; tension and ergonomic adjustments for arms, keyboards & monitors; and integration of new IT or accessory devices.

On each included preventive maintenance visit, our trained techs...

- ✓ Check monitor settings & work surface
- ✓ Check buttons on keyboard for wear and correct operation
- ✓ Test height adjustment operation by raising and lowering work surface
- ✓ Adjust the arm's resistance to counter balance the arm properly
- ✓ Inspect and correct cable placement and management
- ✓ Inspect and tighten bolts, screws, clamps and casters, as applicable

Contact US

(1) 800 947 3901 or repairs@touchpointmed.com



See more before & after photos from recent TPM maintenance visits!

TouchPoint
MEDICAL

20250507-FY-AccessMount Preventive Maintenance-EN